

# FOUNDRY-MASTER Smart 2

## Application Note – Calibration Ranges

# HITACHI

Precision made practical. Confidence in every spark.

### Background

The FOUNDRY-MASTER Smart 2 is a stationary spark spectrometer designed to deliver precise, repeatable, and stable results for routine metal analysis in foundries and manufacturing environments. It is the ideal solution for the metal production and processing industry, providing cost-effective, reliable analysis. This new generation optical emission spectrometer offers high analytical performance at an extremely compact size for seamless quality control in metal production and processing at multiple stages of the process.

Upgraded with CMOS sensor technology, an ultra-stable spark source and improved calibrations covering ten material bases, the FOUNDRY-MASTER Smart 2 ensures excellent repeatability and extended application coverage across ferrous and non-ferrous alloys.

The compact housing and robust construction allow continuous use in demanding industrial conditions.

With the intuitive operating software SpArcfire, straightforward calibration routines, and low maintenance requirements, the FOUNDRY-MASTER Smart 2 provides dependable performance for QA/QC, PMI, and secondary metal production – combining high analytical quality with proven ease of use.

### Key Features

- Analyses the majority of metals and their alloys
- High analytical performance from powerful optics
- Scientific CMOS sensor technology for improved precision
- Ultra-stable spark source for consistent long-term stability
- Wavelength range: 172 – 420 nm
- Very short start-up and measurement time
- Identification of duplex steels with nitrogen analysis
- Compact footprint of just 415 × 665 mm and 35 kg – easy to handle
- Robust housing and design for continuous use in industrial environments
- Preinstalled GRADE Database for fast and easy grade identification
- Results at your fingertips: wide range of result forms, automatic storage
- Excellent price-performance ratio



### \*\*Disclaimer\*\*

The calibration ranges and analytical parameters shown in this document reflect the current development status of the FOUNDRY-MASTER Smart 2. They are not final and may be subject to change as validation and further optimization are completed. Some details regarding element coverage, limits of detection, and precision values may therefore differ in the final release.

# FOUNDRY-MASTER Smart 2

Application Note – Calibration Ranges

# HITACHI

## Sample preparation

Correct sample preparation is very important for precise and accurate OES results. A flat sample surface is essential. To achieve this, different techniques, like grinding or milling, are appropriate, depending on the material and the elements to be analysed.

### Sample preparation for Al, Cu, Mg, Pb, Sn, Ti, Zn

Our recommendation is to use a milling machine equipped with indexable inserts specified for copper alloys. The machine should be optimised for each Cu alloy.

Alternatively, you can use a turning lathe.

For the results presented in this application note, all copper alloys were milled.

### Sample preparation for Co, Fe, Ni

If low Al concentrations need to be measured, aluminum oxide is typically used, depending on the analyte material. Zirconium oxide or silicon carbide are alternatives, grain size 40 – 80.

Cast iron samples are typically prepared with grindstones or cup wheels (stone with segments) while steel is typically prepared with disc or belt grinding machines.

In this case, in order to perform sets of precession measurements, all samples were carefully and appropriately ground on a stationary disk grinder with mesh size 60 Al-corundum paper.



# FOUNDRY-MASTER Smart 2

Application Note – Calibration Ranges

# HITACHI

## Difference between LOD and LOQ


The BEC (equivalent concentration of spectral background) value is the concentration of the analysis sample required to produce the same intensity signal as the background at a given wavelength. The BEC is obtained from the calibration curve and is a fundamental process variable as it directly affects the LOD (**limit of detection**). The LOD is the smallest amount of an element detectable, and it is calculated as follows:

$$LOD = \frac{3}{100} RSD_0 \times BEC$$

$RSD_0$  is correlated to the relative standard value of spectral background. With the BEC value calculated from the calibration curve, we are able to detect different elements in the matrix down to the level of precision (1σ).

However, the **lowest quantitatively determinable amount** (Limit of Quantitation or LOQ) must be larger than the spectrometric LOD by a multiple of three. The resulting LOQ is the quantitatively readable value with our instrument.

The following tables show the calibration ranges of the FOUNDRY-MASTER Smart 2



# FOUNDRY-MASTER Smart 2

Sub-programs & calibration range Fe base

Min and Max in wt%

	Fe_000		Fe_100		Fe_200		Fe_250		Fe_300		Fe_400		Fe_500	
	Orientation		Low alloy		Cast steel		Cr hard & Ni resist		Stainless steel		Tool steel		Mn steel	
	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max
Al	0.0030	6	0.0010	1.5	0.0010	0.28	0.0010	0.28	0.0020	6	0.0015	0.01	0.0020	0.3
As			0.0010	0.2	0.0010	0.2								
B	0.0005	1.8	0.0005	0.07	0.0005	0.1	0.0010	0.1	0.0005	0.015				
Bi			0.0050	0.12	0.0030	0.015								
C	0.0025	4	0.0020	1.6	1.5000	5	1.0000	3.8	0.0030	1.5	0.0060	2.75	0.0040	1.8
Ca			0.0001	0.015					0.0005	0.015				
Ce					0.0050	0.1								
Co	0.0010	10	0.0010	0.55	0.0020	0.25	0.0010	0.6	0.0010	12.5	0.0010	13	0.0010	0.6
Cr	0.0030	35	0.0005	5.5	0.0020	3	0.0015	37.5	0.0020	35.	0.0015	21.5	0.0010	4.2
Cu	0.0020	9	0.0005	1	0.0010	3	0.0010	10	0.0010	8.5	0.0010	0.25	0.0010	0.8
La					0.0005	0.03								
Mg	0.0005	0.1			0.0002	0.2	0.0005	0.2						
Mn	0.0020	20	0.0010	2.5	0.0010	2.5	0.0010	2	0.0010	13.5	0.0015	2	5.0000	25
Mo	0.0050	10	0.0010	1.5	0.0020	1.7	0.0020	4	0.0010	7	0.0025	10.5	0.0015	2.2
N									0.0500	1.2				
Nb	0.0025	2	0.0010	1.6	0.0010	0.17	0.0010	0.75	0.0040	3.3	0.0015	1.7	0.0010	0.35
Ni	0.0030	55	0.0030	5.5	0.0010	5	0.0050	28.5	0.0020	55	0.0010	3.5	0.0020	4.5
P			0.0025	0.12	0.0005	1.5	0.0010	0.73	0.0010	0.15	0.0005	0.35	0.0020	0.12
Pb			0.0030	0.15	0.0040	0.13			0.0080	0.3				
S			0.0010	0.1	0.0020	0.3	0.0020	0.2	0.0020	0.1	0.0020	0.09	0.0020	0.3
Sb			0.0025	0.15	0.0020	0.27			0.0030	0.15				
Se			0.0025	0.4	0.0025	0.03			0.0030	0.4				
Si	0.0020	10	0.0010	2.4	0.0010	6.7	0.0020	6.5	0.0010	4	0.0010	1.3	0.0020	2.5
Sn			0.0010	0.25	0.0010	0.22	0.0010	0.25	0.0010	0.3	0.0010	0.32	0.0010	0.075
Ti	0.0020	2	0.0005	0.5	0.0010	0.26	0.0010	0.27	0.0010	2.5	0.0005	0.5	0.0010	1.2
V	0.0020	10	0.0010	0.8	0.0010	0.75	0.0010	0.75	0.0010	1.	0.0010	10.5	0.0010	0.15
W	0.0050	20	0.0010	1.8	0.0050	0.2	0.0050	0.6	0.0070	3.5	0.0020	22.5		
Zn			0.0005	0.025	0.0010	0.044								
Zr			0.0015	0.4	0.0015	0.05			0.0015	0.4				

# FOUNDRY-MASTER Smart 2

Sub-programs & calibration range Al base

Min and Max in wt%

	AL_000		AL_100		AL_200		AL_300		AL_400		AL_450		AL_500		AL_550	
	Al-orientation		Al-low alloy		Al-Cu alloy		Al-Mg alloy		Al-Si alloy		Al-Si-Cu alloy		Al-Zn alloy		Al-Zn-Si alloy	
	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max
Ag	0.005	1.2	0.0005	1.2	0.0005	1.2	0.0005	1.2	0.0005	1.2	0.0005	1.2	0.0005	1.2	0.0005	1.2
B			0.0005	0.025	0.0005	0.025	0.0005	0.025	0.0005	0.025	0.0005	0.025	0.0005	0.025	0.0005	0.025
Be			0.0001	0.02	0.0001	0.02	0.0001	0.02	0.0001	0.02	0.0001	0.02	0.0001	0.02	0.0001	0.02
Bi			0.0030	0.7	0.0030	0.7	0.0030	0.7	0.0030	0.7	0.0030	0.7	0.0030	0.7	0.0030	0.7
Ca			0.0001	0.02	0.0001	0.02	0.0001	0.025	0.0001	0.02	0.0001	0.02	0.0001	0.02	0.0001	0.02
Cd			0.0005	0.34	0.0005	0.34	0.0005	0.34	0.0005	0.34	0.0005	0.34	0.0005	0.34	0.0005	0.34
Ce			0.0010	0.06	0.0010	0.06	0.0010	0.1	0.0010	0.06	0.0010	0.06	0.0010	0.06	0.0010	0.06
Co	0.005	0.5	0.0010	0.5	0.0015	0.5	0.0010	0.5	0.0010	0.5	0.0010	0.5	0.0010	0.5	0.0010	0.5
Cr	0.005	0.7	0.0005	0.6	0.0010	0.6	0.0005	0.7	0.0005	0.6	0.0005	0.55	0.0005	0.7	0.0005	0.7
Cu	0.010	10	0.0020	0.8	0.0015	10	0.0020	5.5	0.0020	4	0.0020	10	0.0020	4	0.0020	4
Fe	0.005	2.8	0.0015	2.2	0.0010	2.8	0.0015	2.2	0.0015	2.2	0.0015	2.2	0.0015	2.2	0.0015	2.2
Ga			0.0005	0.12	0.0005	0.12	0.0005	0.15	0.0005	0.12	0.0005	0.15	0.0005	0.15	0.0005	0.15
Hg			0.0020	0.075	0.0100	0.075	0.0100	0.08	0.0020	0.075	0.0020	0.075	0.0020	0.075	0.0020	0.075
In			0.0010	0.12	0.0010	0.12	0.0010	0.12	0.0010	0.12	0.0010	0.12	0.0010	0.12	0.0010	0.12
La			0.0005	0.04	0.0005	0.1	0.0005	0.1	0.0005	0.05	0.0005	0.04	0.0005	0.04	0.0005	0.04
Li			0.0002	0.025	0.0002	0.025	0.0002	0.025	0.0002	0.025	0.0002	0.025	0.0002	0.025	0.0002	0.025
Mg	0.005	12	0.0010	1.5	0.0005	1.5	0.0001	12	0.0005	3	0.0010	3	0.0010	4	0.0010	4
Mn	0.005	2.2	0.0005	2.2	0.0005	2.2	0.0005	2.2	0.0005	2.2	0.0005	2.2	0.0015	2.2	0.0015	2.2
Mo	0.005	1	0.0005	1	0.0010	1	0.0005	1	0.0005	1	0.0005	1	0.0005	1	0.0005	1
Na			0.0002	0.015	0.0002	0.015	0.0002	0.015	0.0002	0.015	0.0002	0.015	0.0002	0.015	0.0002	0.015
Ni	0.005	6	0.0010	2.2	0.0010	5.5	0.0010	5.5	0.0010	6	0.0010	5.5	0.0010	6	0.0010	6
P			0.0035	0.018	0.0040	0.018	0.0025	0.018	0.0030	0.018	0.0030	0.018	0.0035	0.018	0.0035	0.018
Pb			0.0010	1.6	0.0010	1.6	0.0010	1.8	0.0010	1.6	0.0010	1.6	0.0010	1.6	0.0010	1.6
Sb			0.0050	0.6	0.0080	0.6	0.0063	0.6	0.0075	0.6	0.0075	0.6	0.0050	0.6	0.0050	0.6
Sc			0.0005	0.85	0.0005	0.85	0.0010	0.85	0.0005	0.85	0.0005	0.85	0.0005	0.9	0.0005	0.9
Si	0.005	25	0.0015	1.5	0.0020	1.7	0.0010	1.35	0.0010	25	0.0015	25	0.0015	1.8	0.0015	25
Sn	0.005	4	0.0020	1.2	0.0020	4	0.0050	4	0.0035	4	0.0020	4	0.0020	4	0.0020	4
Sr			0.0001	0.15	0.0001	0.16	0.0001	0.15	0.0001	0.15	0.0001	0.15	0.0001	0.15	0.0001	0.15
Ti			0.0005	0.65	0.0005	0.65	0.0005	0.65	0.0005	0.65	0.0005	0.65	0.0005	0.7	0.0005	0.7
V			0.0010	0.15	0.0010	0.12	0.0010	0.65	0.0010	0.12	0.0010	0.12	0.0010	0.12	0.0010	0.12
Zn	0.005	12.5	0.0015	2.2	0.0020	2.2	0.0015	5	0.0020	2.2	0.0015	2	0.0020	12.5	0.0020	12.5
Zr			0.0005	1	0.0005	1	0.0005	1	0.0005	1	0.0005	1	0.0005	1	0.0005	1

# FOUNDRY-MASTER Smart 2

Sub-programs & calibration range Cu base

Min and Max in wt%

	Cu_000		Cu_100		Cu_200		Cu_300		Cu_350		Cu_400		Cu_450		Cu_500		Cu_600	
	Global		low alloy copper		Cu-Zn brass		Cu-Sn-Pb bronze		Gunmetal		Cu-Ni alloy		Cu-Ni-Zn alloy		Cu-Al alloy		Be/Co/Ag alloy	
	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max
Ag	0.0015	1.1	0.0015	0.5	0.0015	0.5	0.0015	1	0.0015	1			0.0015	0.12			0.0015	1
Al	0.0020	12	0.0010	1	0.0020	8	0.0010	4	0.0010	0.5	0.0010	1.6	0.0015	0.3	2.0000	12	0.0010	0.2
As			0.0010	0.5	0.0010	0.5	0.0010	1.8	0.0010	0.55					0.0015	0.3		
B			0.0015	0.04								0.0015	0.04					
Be	0.0010	2	0.0002	0.6								0.0001	0.18				0.0002	2.2
Bi			0.0015	0.5	0.0015	6	0.0020	6	0.0015	1	0.0010	0.4						
Cd			0.0005	1	0.0005	0.6	0.0015	0.2	0.0015	0.5								
Co	0.0020	3.3	0.0020	1	0.0020	0.5	0.0015	0.5	0.0020	0.25	0.0015	0.2	0.0015	0.25			0.0010	3.3
Cr	0.0010	3	0.0010	1	0.0010	0.1	0.0015	0.1	0.0015	0.1	0.0010	3					0.0010	0.1
Fe	0.0030	6.7	0.0020	1	0.0020	2	0.0020	2.4	0.0020	1	0.0015	3	0.0020	0.5	0.1000	6.7	0.0020	0.2
Mg			0.0005	0.2							0.0001	0.1			0.0001	0.2		
Mn			0.0005	2	0.0015	20	0.0020	4	0.0015	0.6	0.0015	4.5	0.0015	0.7	0.0015	7	0.0015	0.2
Nb												0.0020	0.8					
Ni	0.0030	36	0.0015	1	0.0015	5	0.0015	3	0.0015	2.2	5.0000	36	5.0000	20	0.0015	8	0.0025	0.5
P			0.0015	0.2	0.0020	0.5	0.0010	1.0	0.0015	0.3	0.0015	0.1	0.0015	0.15			0.0015	0.1
Pb	0.0015	22	0.0015	1	0.0020	5	0.0015	23	0.0015	10	0.0010	0.5	0.0015	1.2	0.0015	0.5	0.0015	0.25
S			0.0010	0.2	0.0010	0.2	0.0010	0.2	0.0010	0.2	0.0015	0.15	0.0015	0.1	0.0010	0.2		
Sb			0.0020	0.5	0.0050	0.5	0.0020	4.5	0.0030	1.2								
Se			0.0010	0.1	0.0015	0.5	0.0020	1										
Si	0.0020	7.8	0.0020	2.2	0.0020	4	0.0020	7.5	0.0020	0.5	0.0015	2	0.0015	1	0.0020	1	0.0020	1
Sn	0.0015	18.5	0.0015	1	0.0015	10	0.0015	18.5	0.0015	11	0.0015	0.5	0.0015	0.5	0.0015	1	0.0015	0.2
Te			0.0040	0.5														
Ti												0.0010	0.075					
Zn	0.0040	50	0.0020	1	0.0030	50	0.0050	12	0.0020	14	0.0030	1.5	10.0000	40	0.0040	1	0.0020	0.25
Zr			0.0005	0.1								0.0010	0.1					

# FOUNDRY-MASTER Smart 2

Sub-programs & calibration range Ni base

Min and Max in wt%

	Ni_000 Global		Ni_100 Low Alloy		Ni_200 Monel		Ni_300 Nimonic/Waspalloy		Ni_400 Incoloy		Ni_500 Inconel		Ni_600 Hastelloy		Ni_700 MarMalloy	
	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max	Min	Max
Al	0.0005	8.5	0.0020	1	0.0010	4.5	0.0030	8	0.0003	0.5	0.0010	8	0.0005	8	1.0000	8
B			0.0001	0.03			0.0010	0.25	0.0003	0.5	0.0003	0.5	0.0005	0.5	0.0005	0.5
C			0.0050	0.5	0.0050	0.5	0.0020	0.5	0.0020	0.5	0.0010	0.5	0.0100	0.5	0.0100	0.5
Co	0.0100	30	0.0010	1	0.0030	2	0.0005	30	0.0005	5	0.0010	30	0.0005	30	2.0000	15
Cr	0.0010	35	0.0020	0.9	0.0020	1	10.0000	24	15.0000	25	10.0000	30	0.0005	25	4.0000	30
Cu	0.0020	38	0.0002	1	22.0000	38.	0.0005	1	0.0010	4	0.0005	3	0.0010	2	0.0010	2.5
Fe	0.0005	55	0.0010	1	0.0005	3	0.0005	35	1.0000	47	0.0005	40	0.0010	25	0.0010	5
Hf	0.0030	2.3													0.0050	2.3
Mg			0.0001	0.15	0.0005	0.25	0.0001	0.2	0.0010	0.25	0.0001	0.25	0.0001	0.2	0.0001	0.25
Mn	0.0100	3.3	0.0010	1	0.0005	3.3	0.0050	1	0.0100	2	0.0050	2	0.0100	2	0.0005	1
Mo	0.0010	40					0.0050	20	0.0500	10	0.0010	12	0.0050	35.	0.0006	4
Nb	0.0100	8.5	0.0030	1			0.0010	8			0.0050	8	0.0100	5	0.0050	2
P			0.0002	0.065	0.0020	0.065										
Pb					0.0020	0.1										
S			0.0002	0.15	0.0005	0.15	0.0010	0.15	0.0011	0.15	0.0010	0.15	0.0010	0.15	0.0010	0.15
Si	0.0005	7.8	0.0010	0.6	0.0010	7	0.0010	3	0.0010	7.8	0.0005	1.5	0.0100	2	0.0010	8
Sn							0.0010	0.8	0.0005	1	0.0025	0.8	0.0010	0.8	0.0020	0.8
Ta	0.0100	7.8													0.0100	7.8
Ti	0.0010	6.5	0.0010	1	0.0005	2	0.0005	6	0.0010	6	0.0005	6	0.0010	5	0.0010	5
V					0.0020	0.5	0.0010	1.2			0.0010	1.2	0.0010	1	0.0050	1.2
W	0.0300	13.5					0.0050	11.5			0.0150	5	0.0150	7	0.0100	13.5
Zr							0.0010	0.3			0.0010	0.3	0.0005	0.3	0.0010	0.3

Min and Max in wt%

	<b>Co_000</b>	
	<b>Global</b>	
	<b>Min</b>	<b>Max</b>
<b>Al</b>	0.003	1.6
<b>B</b>	0.003	0.2
<b>C</b>	0.003	2.8
<b>Cr</b>	5.000	36
<b>Cu</b>	0.002	0.2
<b>Fe</b>	0.005	21
<b>Mn</b>	0.005	2.1
<b>Mo</b>	0.003	8.7
<b>Nb</b>	0.005	2.8
<b>Ni</b>	0.010	26
<b>P</b>	0.003	0.05
<b>Pb</b>	0.006	0.045
<b>S</b>	0.003	0.07
<b>Si</b>	0.005	1.8
<b>Sn</b>	0.002	0.15
<b>Ta</b>	0.005	4.2
<b>Ti</b>	0.002	0.5
<b>W</b>	0.005	17
<b>Zr</b>	0.003	0.5

## FOUNDRY-MASTER Smart 2

Sub-programs & calibration range Co base

Min and Max in wt%

	<b>Mg_000</b>	
	<b>Global</b>	
	<b>Min</b>	<b>Max</b>
<b>Ag</b>	0.0010	3
<b>Al</b>	0.0050	13.8
<b>Be</b>	0.0001	0.004
<b>Ca</b>	0.0005	0.7
<b>Cd</b>	0.0005	0.08
<b>Ce</b>	0.0050	2.3
<b>Cu</b>	0.0020	3.2
<b>Er</b>	0.0005	0.25
<b>Fe</b>	0.0050	0.04
<b>Gd</b>	0.0030	2.3
<b>La</b>	0.0020	1.5
<b>Mn</b>	0.0005	2.5
<b>Nd</b>	0.0090	3.5
<b>Ni</b>	0.0025	0.06
<b>Pb</b>	0.0050	0.1
<b>Pr</b>	0.0050	0.9
<b>Si</b>	0.0200	2
<b>Sn</b>	0.0020	0.1
<b>Sr</b>	0.0005	1.2
<b>Yb</b>	0.0010	0.15
<b>Zn</b>	0.0100	7
<b>Zr</b>	0.0010	0.6

## FOUNDRY-MASTER Smart 2

Sub-programs & calibration range Mg base

Min and Max in wt%

	<b>Ti_000</b>	
	<b>Global</b>	
	<b>Min</b>	<b>Max</b>
<b>Al</b>	0.005	9.00
<b>C</b>	0.001	0.13
<b>Cr</b>	0.001	7.00
<b>Cu</b>	0.002	3.00
<b>Fe</b>	0.005	2.50
<b>Mn</b>	0.002	7.70
<b>Mo</b>	0.010	5.30
<b>Nb</b>	0.010	7.50
<b>Ni</b>	0.001	1.00
<b>Pd</b>	0.001	0.20
<b>Ru</b>	0.005	0.11
<b>Si</b>	0.002	0.70
<b>Sn</b>	0.005	12.40
<b>Ta</b>	0.005	1.10
<b>V</b>	0.005	16.20
<b>W</b>	0.010	1.10
<b>Zr</b>	0.005	6.00

## FOUNDRY-MASTER Smart 2

Sub-programs & calibration range Ti base

# FOUNDRY-MASTER Smart 2

Sub-programs & calibration range Zn base

Min and Max in wt%

	<b>Zn_000</b>		<b>Zn_100</b>		<b>Zn_200</b>		<b>Zn_300</b>		<b>Zn_400</b>		<b>Zn_500</b>	
	<b>Orientation</b>		<b>low alloy</b>		<b>Zn-Al 2-6%</b>		<b>Zn-Al 6-15%</b>		<b>Zn-Al &gt;15%</b>		<b>Zamak</b>	
	<b>Min</b>	<b>Max</b>	<b>Min</b>	<b>Max</b>	<b>Min</b>	<b>Max</b>	<b>Min</b>	<b>Max</b>	<b>Min</b>	<b>Max</b>	<b>Min</b>	<b>Max</b>
<b>Al</b>	0.0010	60	0.0005	1.5	1.5000	6.2	5.5000	15.5	15.0000	60	3.0000	5
<b>Bi</b>			0.0005	0.15								
<b>Cd</b>	0.0005	0.8	0.0005	0.8	0.0005	0.04	0.0005	0.06	0.0005	0.1	0.0005	0.03
<b>Cr</b>			0.0010	0.2								
<b>Cu</b>	0.0010	6.5	0.0005	1	0.0010	3.5	0.0015	6.5	0.0020	3.5	0.0005	4
<b>Fe</b>	0.0050	0.7	0.0020	0.7	0.0015	0.12	0.0015	0.5	0.0015	0.7	0.0015	0.03
<b>Mg</b>	0.0001	0.25			0.0001	0.25	0.0001	0.15	0.0001	0.15	0.0001	0.1
<b>Mn</b>	0.0010	0.1			0.0010	0.1	0.0010	0.02	0.0005	0.015		
<b>Ni</b>					0.0020	0.05	0.0020	0.015			0.0015	0.03
<b>Pb</b>	0.0010	2.5	0.0005	3	0.0010	0.05	0.0005	0.08	0.0010	2.5	0.0005	0.03
<b>Sb</b>			0.0025	0.4								
<b>Si</b>	0.0010	0.07			0.0010	0.05	0.0005	0.07	0.0005	1.2	0.0005	0.03
<b>Sn</b>	0.0010	3	0.0005	3	0.0015	0.2	0.0010	0.15	0.0010	0.02	0.0010	0.015
<b>Ti</b>							0.0005	0.35				

Min and Max in wt%

	<b>Pb_000</b>	
	<b>Global</b>	
	<b>Min</b>	<b>Max</b>
<b>Ag</b>	0.0005	2.5
<b>Al</b>	0.0005	0.6
<b>As</b>	0.0010	1.6
<b>Au</b>	0.0005	0.1
<b>Ca</b>	0.0005	0.27
<b>Cd</b>	0.0010	0.5
<b>Cu</b>	0.0005	0.4
<b>Fe</b>	0.0010	0.02
<b>In</b>	0.0005	0.7
<b>Ni</b>	0.0005	0.02
<b>Sb</b>	0.0010	11
<b>Se</b>	0.0010	0.05
<b>Sn</b>	0.0010	70
<b>Te</b>	0.0010	0.04
<b>Tl</b>	0.0005	0.06
<b>Zn</b>	0.0005	0.5

## FOUNDRY-MASTER Smart 2

Sub-programs & calibration range Pb base

Min and Max in wt%

	<b>Sn_000</b>		<b>Sn_200</b>	
	<b>Global</b>		<b>lead free solder</b>	
	<b>Min</b>	<b>Max</b>	<b>Min</b>	<b>Max</b>
<b>Ag</b>	0.0005	4.7	0.0005	4.7
<b>Al</b>	0.0002	0.075	0.0002	0.075
<b>As</b>	0.0050	0.6	0.0050	0.6
<b>Au</b>	0.0005	0.125	0.0005	0.125
<b>Bi</b>	0.0005	55	0.0005	55
<b>Cd</b>	0.0005	2	0.0005	0.55
<b>Co</b>	0.0005	0.025	0.0005	0.025
<b>Cu</b>	0.0005	12	0.0005	12
<b>Fe</b>	0.0005	0.08	0.0005	0.08
<b>Ga</b>	0.0005	0.05	0.0005	0.05
<b>Ge</b>	0.0050	0.5		
<b>Hg</b>	0.0010	0.17		
<b>In</b>	0.0005	0.15	0.0005	0.15
<b>Ni</b>	0.0005	1.4	0.0005	1.4
<b>P</b>	0.0020	0.035	0.0020	0.035
<b>Pb</b>	0.0005	40	0.0005	1
<b>S</b>	0.0020	0.025	0.0020	0.025
<b>Sb</b>	0.0005	15.5	0.0005	15.5
<b>Se</b>	0.0005	0.05		
<b>Te</b>	0.0005	0.08		
<b>Zn</b>	0.0005	9	0.0005	9

## FOUNDRY-MASTER Smart 2

Sub-programs & calibration range Sn base

# FOUNDRY-MASTER Smart 2

Application Note – Calibration Ranges

# HITACHI

## PERFORMANCE DISCLAIMER

All performance data provided in this document is based on certified reference materials and typical sample conditions (flat ground or milled). The calibration ranges shown represent averaged data across multiple material types and should be considered as indicative "typical" values only. Calibration ranges can be extended upon request using verified customer samples.

Detection limits (DL) have been determined on a production series instrument according to **DIN 32645:2008-11** using a 99.7% confidence interval ( $3\sigma$  method), and reflect ideal conditions with well-prepared, homogeneous reference materials. Please note that these limits may vary depending on the alloy matrix, sample homogeneity, and measurement parameters. **A reported value at or near the detection limit should not be interpreted as definitive evidence of the presence of an element unless the matrix and signal behavior clearly support it.**

Upper and lower concentration ranges are influenced by the availability and suitability of certified or secondary standards. Extensions beyond published values are possible if sufficient suitable standards are available and verified for homogeneity and traceability.

For specific applications or material requirements, we recommend discussing your needs directly with our technical team. We are happy to verify performance with your own samples upon request.

# FOUNDRY-MASTER Smart 2

Application Note – Calibration Ranges

# HITACHI

## Summary

Hitachi High-Tech's FOUNDRY-MASTER Smart 2 delivers accurate and repeatable results for a wide range of ferrous and non-ferrous alloys. Its compact footprint, robust design and intuitive operation make it well suited for routine use in foundries and manufacturing environments. The combination of CMOS sensor technology, the proven spark source, and extended calibrations ensures stable performance across multiple applications, from QA/QC to secondary metal production. With results available within seconds and automatic data storage, the FOUNDRY-MASTER Smart 2 supports fast decision-making, improved process control, and reduced production costs.

## Ordering information

The instrument configuration and accessories used to produce the data in this application note are included in the following packages:

FOUNDRY-MASTER Smart 2 with standard calibrations (10 bases) or customer-specific calibration set

### Optional accessories

- Wire adapter kit (1.5–10 mm diameter)
- Boron nitride sample plate for small parts ( $\geq 5$  mm)
- Additional sample holders and spark stand inserts

## More information

To find out more visit: [www.hitachi-hightech.com/hha](http://www.hitachi-hightech.com/hha)

This publication is the copyright of Hitachi High-Tech Analysis Corporation and provides outline information only, which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or regarded as the representation relating to the products or services concerned. Hitachi Hitachi High-Tech Analysis Corporation's policy is one of continued improvement. The company reserves the right to alter, without notice the specification, design or conditions of supply of any product or service. Hitachi High-Tech Analysis Corporation acknowledges all trademarks and registrations.

© Hitachi High-Tech Analysis Corporation, 2025 All rights reserved.

